

### About Us:

The Goal of GetMeKnown is to create an Online Microchip Management Service that will allow you to manage all your registered Micro-Chips in an easy to understand manner while still giving you complete control.

With features that allow you to transfer ownership (no need to waste time to contact us), complete with logging of all actions as evidence/history, etc.

You're in Good Hands and here's why;

- x Get Me Known adheres to the International standards like ISO11784/5 FDX-B, ICAR, etc. Read by all ISO animal scanners.
- x Backed by a Company in existence and supporting "dogdom" for the past 10 years with various other products.
- x International Accessible Database stored on rack servers - We have never been offline from public launch from inception.
- x Innovation of a breeder to assist the breeders, shelters, Vets, for Import and Export, etc. No hidden cost.
- x Extra features like to flag a missing animal, Record card for shelter management, etc..
- x **No Annual Fees !!!!!!!!**

### Our Goal :

To make micro chipping affordable with high quality products to allow each responsible owner, Vet or Shelter to microchip and register their pets/animals online and "get them known" when lost or legally to verify identity!!

### Register a New Account

To perform any account related action on **GetMeKnown** (<http://www.getmeknown.co.za/>) you'll need an account, and creating one is quite easy. To start head over to <http://www.getmeknown.co.za/account/register> and start the process. Reminder that you are the injector and to use your id no. or your KUSA/Other membership number.

### Register a Microchip

After you have created a account you will need to Login (<http://www.getmeknown.co.za/account/login>), this is done by providing the E-Mail Address you used to register and the password for your account. To create a Microchip click on "**Register Microchip**" and follow the directions on the next screen. You will then be asked to Provide your Microchip Number. A Microchip Numbers consists of 15 Digits. Please verify correctness. Complete all the required fields. View all your registered pets –just click "account".

### Print or e-mail Certificates

Once you have completed a registration of a microchip you can view, print or e-mail the certificate. Should you need a copy at a later stage just login to your account and press "Certificate" next to the desired microchip number.

### Transfer ownership of a microchip

To Receive Ownership of a Microchip you will need to have an account and Request Ownership and wait for the current owner to login and then approve your request. **No need to type in any information about the Pet.** Login to your account, then in the Main Menu on your profile page select "Request Ownership". Type in the 15 digit microchip number. You will be redirected to a page where you can review the details of the pet click "Request a Transfer of Ownership" the request will be sent you just need to wait for the current owner to approve.

## RFID MICROCHIPS FOR DOGS, CATS, CATTLE, SHEEP, WILDLIFE, ETC...


Identify animals with an unique code, to prove ownership, stud registrars, etc. ICAR – International Recorded Numbering & ISO protocols 11784/85 FDX-B. Read by all ISO Veterinary, Wildlife and Rescue organization scanners (15digit). Each Microchip comes with a sterile disposable Injector and Barcode Stickers

Visit [www.getmeknown.co.za](http://www.getmeknown.co.za) or send us a e-mail at [info@lusahn.co.za](mailto:info@lusahn.co.za) Tel: 0823736375 (Adri )

# GetMeKnown (PTY) Ltd 2015/436468/07 t/a [www.GetMeKnown.co.za](http://www.GetMeKnown.co.za)

O Box 1081, Malmesbury, 7299 [info@lusahn.co.za](mailto:info@lusahn.co.za) Tel 0224857310 / Fax 088 0224857310

“Reputation for Service Excellence”

<p><b>RFID (microchip) Concept:</b> This is not a GPS microchip and cannot be tracked via satellite. It is a permanent identification chip that is implanted just under the animals’ skin between the shoulder blades in the neck and must be scanned with a microchip reader/scanner at close range (max +-10cm). You must register the microchip online and then when the database is searched you are identified as the owner. Every Animal microchip have an unique identification..</p> 	<p><b>Our Products:</b> <b>Microchips:</b> Our Microchips adhere to the ISO Standard 11784/85 FDX-B 15 digit international standards and is also ICAR recorded (International 15 digit Microchips). Our Microchips is only 1.4mmx8mm. Smallest around.</p> <p><b>Readers:</b> The Lusahn Halo RFID Reader are able to read all the 15 digit (FDX-B), older 10 digit ID codes (FDX-A) and even the Bio-Thermo Microchips. Imported from the UK where its extensively used by Vets, Police, shelters, breeders, etc. No excuses anymore we read them ALL ...</p>		
<p><b>Our Database:</b> International Accessible database on rack servers <b>We have never been offline from public launch.</b> Our Infrastructure is Maintained by Certified Engineers. No paperwork as you manage all your microchips online, you can even e-mail certificate to all relevant people. Accepted by KUSA, CASA, pet medical aids, etc. You can even flag a lost animal to display extra info on the screen</p>	<p><b>Contact Us:</b> Website: <a href="http://www.getmeknown.co.za">www.getmeknown.co.za</a> Office: 022 485 7310 Cell: 082 373 6375 Fax: 0880224857310 fax to email e-mail: <a href="mailto:info@lusahn.co.za">info@lusahn.co.za</a> <b>Contact : Adri/Johan</b></p>		
<b>Microchip Prices</b>			
1 – 4	R70.00	20 – 49	R55.00
5 – 9	R65.00	50 – 99	R50.00
10 – 19	R60.00	100 – more	R45.00
<b>NPO Rescues Only &gt;100 order</b>	<b>R40.00</b>		
Halo Scanner FDX (for all international 15 and 10 Digits Microchips/Transponders)			R1780

## Unique features :

- You as account holder/ owner may flag the microchip to display a online message when dog was stolen, lost, etc..
- No Annual fees or any extra hidden cost. Life Long Warranty.

## MICROCHIPS FOR DOGS, CATS, CATTLE, SHEEP, WILDLIFE, ETC.

Identify animals with an unique code, to prove ownership, stud registrars, etc.  
ICAR – International Recorded Numbering & ISO protocols 11784/85 FDX-B.  
Read by all ISO Veterinary, Wildlife and Rescue organization scanners (15digit).  
Each Microchip comes with a sterile disposable Injector and Barcode Stickers.

> 100 @ R45
Scanner @ R1780
Quality & Affordable

Visit [www.getmeknown.co.za](http://www.getmeknown.co.za) or send us a e-mail at [info@lusahn.co.za](mailto:info@lusahn.co.za) Tel: 0823736375 ( Adri )